



INDIVIDUAL PMMA guide

Drilling templates for guided implantology

BENDING STRENGTH*

> 105 MPa

INDIVIDUAL PMMA guide



Order Number 001-6689EN

INDICATIONS

Production of drilling templates for guided implantology.

VORTEILE

- + Extremely high hardness and tensile strength values
- + Increased safety through guided implantology process

REFERENCES

Provided as an open STL file from design software (e.g. 3Shape, exocad, or similar) or from inLab SW as a i.lab-file.

The drilling templates are free from supports, without guide sleeves and are non-sterile upon delivery. Dentavon® Liquid (manufacturer: Schüller & Mayr), an aldehyde-free disinfectant based on active oxygen, can be used as a disinfectant.

COMPOSITION*

PMMA		100.0 %
(Polymethylmethacrylat)		

* typical material characteristics.
The published information of the respective material manufacturer shall apply.