

INFINIDENT

SOLUTIONS



Lab Model SL/SLS

High-grade working models for the dental lab

Lab Model SL/SLS



Order Number 001-2080EN

INDICATIONS

As the basis for crown and bridge frameworks in the anterior and posterior regions for up to 16 units.

BENEFITS

- + Sturdier than plaster models: does not shatter if dropped
- + More abrasion-resistant than plaster models: stump shape does not change even if the restoration is inserted and removed many times
- + Already segmented: each prepared stump is its own segment
- + Already pinned to the model plate, similar to the conventional model process
- + Optional additional stump with gingival mask available

REFERENCES

Lab Model SL/SLS* can be used to produce physical working models of acrylate plastic based on the digital impression data from Dentsply Sirona Connect. You have the choice between a single-tooth model (SLS) with three teeth per jaw or the full jaw model with up to 16 teeth per jaw (SL). When you order a single-tooth model you will receive three teeth per jaw (prepared stump plus medial and distal adjacent teeth) regardless of how many teeth a digital impression was made of. When you order the full jaw model, the model corresponds to the digital impression provided by the dentist. Both model types consist of the saw-cut model and the opposite jaw model pinned on a base plate.

Note: The terminal occlusion holder is not included.

* typical material characteristics.
The published information of the respective material manufacturer shall apply.